

Application News

Food - MOSH/MOAH

Mineral Oil Residues in Food Part 1- Dry Food (Flour, Noodles and Rice)

No. SCA-100-001

Introduction

Mineral oil (MO) residues in food raised public concern due to some elevated concentrations up to several thousand milligrams per kilogram food [1]. Due to the chemical structures two groups of MOs can be differentiated. Mineral oil saturated hydrocarbons (MOSH) consist of linear and branched alkanes, and alkylsubstituted cyclo-alkanes, whilst mineral oil aromatic hydrocarbons (MOAH) include mainly alkyl-substituted polyaromatic hydrocarbons. Technical grades of mineral oil contain aromatic hydrocarbons in a concentration range from 15-35%. Food grade mineral oils (white oils) can have lower MOAH concentrations. Rice and pasta and other grain based food products including flour from different grain varieties are consumed in high amounts all over the world. The annual rice consumption per capita is over 54 kg per capita [2]. The annual pasta consumption in Italy is 26 kg and approximately 8 kg in Germany [3].

At the moment there are no legal limits in Europe for MOSH/MOAH, but MOSH concentrations up to 2 mg/kg and MOAH levels below 0,5 mg/kg are considered to be acceptable.

■ System Setup

For the determination of the two mineral oil fractions an online LC-GC –FID system was

used. The LC is directly connected to two high temperature GC columns with retention gaps which are installed in one GC oven. MOSH and MOAH fractions are separated on a silica gel column using a n-hexane /dichloromethane gradient. The interface between LC and GC is controlled by Axel Semrau LC-GC Chronect interface. After transferring the MOSH fraction on column 1 and MOAH on column 2 the temperature programme is started and both fractions are separated simultaneously and detected by FID. Figure 1 shows a typical LC-Chromatogram with UV-signal in black, pump pressure in green, CH₂CI₂ concentration in blue and total flow in purple. Figure 2 shows the LC-GC-FID system.

LC Parameters:

Shimadzu LC-20AD solvent delivery pump Column: Allure Silica 5 μ m (250 × 2.1 mm) Gradient: Start with 100 % n-Hexane (flow 0.3 ml/min), CH₂Cl₂ raised to 35 % within 2 min (hold for 4.20 min), column was flushed at 6.30 min with 100 % CH₂Cl₂ (flow 0.5 ml/min; hold for 9 min) and reconditioned to 100 % n-Hexane (flow 0.5 ml/min; hold for 10 min). Flow was decreased afterwards to 0.3 ml/min until next injection.

UV-Detector: D₂-lamp; 230 nm, 40 °C cell temperature




Fig. 1: LC-Chromatogram

GC Parameters:

Shimadzu GC-2010 Plus dual FID

Guard Columns: Restek MXT Siltek (10 m ×

0.53 mm id)

Columns: Restek MTX ®-1 (15 m × 0.25 mm id

 \times 0.1 μ m df)

Carrier gas: Hydrogen (150 kPa analysis pressure; evaporation pressure: 87 kPA

MOSH, 85 kPA MOAH Temperature program:

60 °C (6 min) @ 20 °C/min to 100 °C (0 min) and followed by 35 °C/min to 370 °C (9.29 min)




Fig. 2: LC-GC-FID System

■ Experimental Work

Sample preparation

Depending on the expected mineral oil concentration 1-10 g of homogenized and finely ground samples were used. The sample was extracted in Hexane/Ethanol 1/1 after the addition of an internal standard mixture (Restek MOSH/MOAH standard Cat.#: 31070 containing 9 internal standards) at room temperature for 2 hours under

occasionally shaking the flask. After the extraction water was added and centrifuged. The washing step was repeated a second time. Afterwards the organic layer was dried over anhydrous sodium sulphate and the volume was adjusted to 1 mL in an automated solvent concentrator. The extract was transferred into a 2 mL autosampler vial and put in the autosampler rack of the LC-GC system. Aliquots of 50 μ l were injected into the LC. Here the separation into the MOSH-and MOAH fraction is performed. Each fraction has a size of 450 μ L and is transferred to the respective channel of the GC for parallel MOSH/MOAH determination.

Quantification

For the quantification of the MOSH fraction the internal standard Bicyclohexyl (CyCy) is used. For the MOAH fraction 2-Methylnaphtahlene (2-MN) is used. All other internal standards are used to ensure no losses of analytes and a good separation between the MOSH and MOAH fraction.

According to a proposed method published by the German Bundesinstitut für Risikobewertung (BfR) quantification is done by integration of the hump for different molecular weight regions. They propose for food contact materials three ranges for the MOSH fraction (C_{10} - C_{16} , C_{16} - C_{25} and C_{25} - C_{35}) and two ranges for the MOAH fraction (C₁₀-C₂₅, C₂₅-C₃₅). For dry food only the ranges up to C₂₅ are used [4]. Figure 3 shows the MOSH (black trace) and MOAH (purple trace) of a spaghetti sample with a MOSH concentration of 12.7 mg/kg from C₁₆-C₃₅ and the marked regions of C₁₆-C₂₅ (blue) and C₂₅-C₃₅ (red). The internal standards are marked with symbols (black squares internal standards MOSH: C₁₁, Bicyclohexyl, C₁₃ and Cholestane eluting with the MOSH hump, purple stars internal standards MOAH: Pentylbenzene, 1 & 2-Methylnaphthalene, Tri-tert-butylbenzene and Perylene) .The rice sample (Figure 3) shows also additional peak in the rear part of the chromatogram. These are naturally occurring odd-numbered alkanes with a chain length of C21 to C35.

Fig. 3: Spaghetti sample




Fig. 4: Rice sample




Fig. 5: n-Alkane standard C10-C40 as retention time marker

■ Conclusion

The new LC-GC-FID provides a tool for fast and reliable routine analysis of MOSH and MOAH in dry non fatty food. After the extraction step the samples are analysed fully automated with a high sample through put.

■ References

4)

- EFSA Panel on Contaminants in the Food Chain (CONTAM) Scientific Opinion on Mineral Oil Hydrocarbons in Food DOI: 10.2903/j.efsa.2012.2704
- 2) Link: https://www.statista.com/statistics/256002/global-per-capita-rice-use-since-2000/
- Link: http://www.internationalpasta.org/index.aspx? idsub=118
 - Link: http://www.bfr.bund.de/cm/343/bestimmungvon-kohlenwasserstoffen-aus-mineraloeloder-kunststoffen.pdf



For Research Use Only. Not for use in diagnostic procedures. Not available in the USA, Canada, and China. This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.