

Application News

No. 277A

Gas Chromatography Mass Spectrometry

High Sensitivity Analysis of Coffee Aroma Components Using the SPME Arrow

Gas chromatograph mass spectrometers (GC-MS) capable of excellent qualitative measurements are used in the analysis of aroma components in foods and beverages. The convenient sampling methods of SPME (solid-phase microextraction) and HS (headspace extraction) are increasingly used for sample introduction. However, sample introduction methods such as these can suffer from insufficient sensitivity when analyzing some aroma components.

The SPME Arrow method was developed as a new sample introduction option for the AOC-6000 Multifunctional Autosampler to address this shortcoming. The larger sorption phase volume compared to conventional SPME fibers allows the SPME Arrow to achieve high enrichment of volatile components and serve as a solution for applications where sensitivity was previously lacking.

This article presents the results of analyzing coffee aroma components using the SPME Arrow.

K. Kawamura

Sample Introduction Using the SPME Arrow

The SPME Arrow enables analysis at high sensitivities due to approximately 5 to 20 times more sorption phase than conventional SPME fibers (Fig. 1).

The thick and robust structure of the SPME Arrow also provides higher durability over conventional SPME fibers.




Fig. 1 Comparison of SPME Arrow and Conventional SPME Fibers

Sample and Analysis Conditions

A sample of 2 g of commercially-available ground coffee beans were weighed out, placed into a 20 mL vial, and set in the AOC-6000. Table 1 lists the analysis conditions. For comparison, analysis was also performed using the conventional SPME method.

Table 1 Analysis Conditions

| | | | |
|--------------------------------|--|------------------------|---|
| GCMS | : GCMS-QP2020 | | |
| Autosampler | : AOC-6000 | | |
| Column | : DB-WAXetr (length: 60 m, 0.25 mm I.D., df = 0.25 µm) | | |
| SPME Arrow conditions | | | |
| SPME Arrow | : PDMS (O.D.: 1.1 mm, film thickness: 100 µm, length: 20 mm) | SPME conditions | |
| Conditioning Temp. | : 270 °C | SPME fiber | : PDMS (film thickness: 100 µm, length: 10 mm) |
| Pre Conditioning Time | : 5 min | Conditioning Temp. | : 270 °C |
| Incubation Temp. | : 60 °C | Pre Conditioning Time | : 5 min |
| Incubation Time | : 8 min | Incubation Temp. | : 60 °C |
| Stirrer Speed | : 250 rpm | Incubation Time | : 8 min |
| Sample Extract Time | : 30 min | Agitator Speed | : 250 rpm |
| Sample Desorb Time | : 2 min (250 °C: GC vaporizing chamber temperature) | Sample Extract Time | : 30 min |
| | | Sample Desorb Time | : 2 min (250 °C: GC vaporizing chamber temperature) |
| GC conditions | | MS conditions | |
| Vaporizing chamber temperature | : 250 °C | Interface temperature | : 230 °C |
| Injection mode | : Splitless | Ion source temperature | : 200 °C |
| Purge flow rate | : 5.0 mL/min | Ionization method | : EI |
| Control mode | : Linear Velocity (25.5 cm/sec) | Measurement mode | : Scan |
| Column oven temperature | : 40 °C (3 min) → 10 °C/min → 250 °C (10 min) | Event time | : 0.3 sec |

Analysis Results

Fig. 2 shows the analysis results for the SPME Arrow and SPME fibers. Known aroma components of coffee that were detected include short-chain aldehydes, phenols, pyridines, and pyrazines in addition to multiple sulfur-based compounds. Compared to the conventional SPME fiber, the SPME Arrow enabled high enrichment and analysis of the aroma components. Fig. 3 shows a comparison of the mass chromatogram peak area of several aroma components.

Conclusion

In contrast to conventional SPME fibers, analysis of trace components that prove difficult with SPME can be achieved using the SPME Arrow, which is coated with a larger volume of sorption phase.




Fig. 2 Comparison of SPME Arrow and Conventional SPME Fibers




Fig. 3 Peak Area Comparison of Typical Aroma Components (Area Comparison Where SPME Arrow Value is 1)